


INTEGRATED WIDE ANGLE VIEWING SYSTEM for FUNDUS OBSERVATION



Features

- Integrated inverter which minimizes the height of the microscope and flips away when not needed
- Ultra view non inverter optical system
- Easy to adapt with any microscope
- Excellent field of view with better stereoscopic vision
- High optical resolution with focusing mechanism to easy the focus of the retina
- Optimum rotation without retinal shift
- Uninterrupted access to surgeon at high diopter






Technical data

- Microscope working distance - 175 mm or 200mm
- Front Mid field lens - for approx. 90° viewing angle
- Front Wide Angle lens - for approx. 124° viewing angle

Distributed by



Mfg by:

AEON MEDITEC (INDIA) PRIVATE LIMITED

Plot No. 13, Door No. 4, 5th Cross Street, 2nd Main Road,

Sri Lakshmi Nagar, Valasaravakkam, Chennai - 600087

www.aeonmeditec.com | info@aeonmeditec.com |

(+91-44) 42840149 | +91-8012531317

Disclaimer: Specifications are subject to change, without prior notice.

The information contained in corresponds to the state of art and alterations.

