



Base Tray

Transparent Lid

Features:

- ❖ Frequency 42kHz / 170W.
- ❖ Heating power 80W.
- ❖ 5 steps digital timer.
- ❖ Built-in heating option.
- ❖ Strong plastic housing: Better water-proof and drip-proof ability than metal.
- ❖ With moulded base tray.
- ❖ With removable transparent lid.

Model	Product	Tank Size (Cap)	App. Shipping Details Net / Gr. / Volume
CD-4820	Digital Ultrasonic Cleaner	24x14x10cm (2.5L)	2.5 kg / 3.5 kg / 0.020 m ³
Optional Accessories			
	Liquid Detergent Concentrate	1 lit bottle	1.0 kg / 1.5 kg / 0.005 m ³

All electrical appliances work on 230V AC, 50Hz, single phase. Other Voltages available on request.
Due to continuous improvement in design, the product supplied may have modified features.