

## Features

- ★ Produces 100% distilled, pure safe water with quality guarantee.
- ★ Output - 4 liters in 4.5 hours.
- ★ Easy to operate, easy to clean.
- ★ Mobile, durable and a space saver.
- ★ Power consumption - 580W
- ★ Dual overheat protection - Thermostatic cut off and Thermal Fuse cut off.
- ★ Power automatically switches off after water distillation is complete.
- ★ Powder coated body with steel inner parts.
- ★ CE approved.
- ★ More than 200 units installed all over India.



IDEAL FOR - Automatic Autoclaves, for cleaning of instruments,  
Ultrasound cleaner and fogging machine

### For Details Contact :

**M** **MEDISYS ENTERPRISES**  
AN ISO 9001: 2008 CERTIFIED CO.

Off No.3, Nair Compound, Film City Road, Goregaon(E), Mumbai - 400062  
Ph.: 91-22-28495501, 28497725 E-mail : admin@medister.in  
Website : www.medisysenterprises.com, www.medister.in

### Our Sales / Service Network :

Maharashtra / New Delhi / West Bengal / Tamilnadu  
Karnataka / Andhra Pradesh / Kerala / Madhya Pradesh  
Gujarat / Rajasthan / Punjab / Uttar Pradesh / Orissa.

