

### Your Partner for Sterilization and Infection Control



# Flash steam sterilizers for rapid sterilization with vacuum and printer



**Mediclave 8B (Flash)**



**Mediclave 22B (Eco)**



**Mediclave 22B (LCD)**

- » **Fast sterilization**
- » **Fully Automatic**
- » **8 Preset Sterilization Programs**
- » **Ster. Temp.: 105-134C**
- » **Ster. Time : 4-60 min.**
- » **Dry Time : 1-25 min.**
- » **Vacuum Times : 1-10 times**
- » **B&D / HELIX / VACUUM Testing**
- » **Cleaning Program**
- » **Printer & USB Output Available**
- » **Big Stainless Steel Water Tanks**
- » **High-quality Bacteria Filter**

## TECHNICAL DATA :

Model	Unit	Mediclave 8B (Flash)	Mediclave 22B (Eco)	Mediclave 22B (LCD)
Capacity	Ltrs.	8	22	22
Inner dimensions	mm	Ø 170 x 320	Ø 250 x 450	Ø 250 x 450
Outer dimensions (WxDxH)	mm	450 x 610 x 430	450 x 650 x 430	500 x 650 x 430
Loading type	-	Front	Front	Front
Heater power	kW	1.8	2.0	2.0
Trays	nos	3	3	3
Special Cycles	-	Vacuum / BD Test / Helix	Vacuum / BD Test / Helix	Vacuum / BD Test / Helix
Operating Temp. / Pressure	°C / kPA	121°C / 134°C	121°C / 134°C	121°C / 134°C
Duty Cycle	-	Eight cycles in an eight hour shift	Eight cycles in an eight hour shift	Eight cycles in an eight hour shift
Min Installation Size	mm	440 x 360	540 x 415	540 x 415
Total Cycle Time	min	14 min (Fastest)	38 min (Fastest)	38 min (Fastest)
Cycle Record	-	USB / Printer	USB / Printer	USB / Printer
Packing Dimensions (WxDxH)	mm	650 x 450 x 420	770 x 570 x 510	770 x 570 x 510
Net / Gross Weight	kg	35 kg	55 kg	55 kg

**For Details Contact :**

## MEDISTER

**Infection Control Systems Pvt. Ltd.**

**Office :** G-119, Gr. Floor, Shagun Arcade Premises, Near Dindoshi Bus Depot, Film City Road, Malad (East), Mumbai - 400 097.  
**Factory :** 312/313, Sarita A, Prabhat Industrial & Res. Complex, Opp. Toll Naka, W. E. Highway, Dahisar (E), Mumbai - 400 068.  
**Phone :** 91-22-28495501 / 28497725, **Email :** admin@medister.in,  
**Website :** www.medister.in



**Our Sales / Service Network :**

Maharashtra / New Delhi / West Bengal / Tamilnadu / Karnataka / Andhra Pradesh / Kerala / Madhya Pradesh / Gujarat / Rajasthan / Punjab / Uttar Pradesh / Orissa / Assam / Himachal / Uttarakhand / Goa / Bihar / Jharkhand / Nagpur / Chattisgarh.