

Your Partner for Sterilization
and Infection Control



Mediclave - Pre and Post Vacuum Autoclaves

Flash steam sterilizers for rapid sterilization with vacuum and printer



Mediclave 12B (Flash)

- » Fast sterilization
- » Fully Automatic
- » 8 Preset Sterilization Programs
- » Ster. Temp.: 105-134C



Mediclave 35B / 50B

- » Ster. Time : 4-60 min.
- » Dry Time : 1-25 min.
- » Vacuum Times : 1-10 times
- » B&D / HELIX / VACUUM Testing



Mediclave 17B / 23B

- » Separate Tank - Clean Water / Used Water
- » Printer & USB Output Available
- » High-quality Bacteria Filter (HEPA)
- » SS Door Cover & Control Panel

TECHNICAL DATA :

Model	Unit	Mediclave 12B (Flash)	Mediclave 17B	Mediclave 23B	Mediclave 35B	Mediclave 50B
Capacity	Ltrs.	12	17	23	35	50
Inner dimensions	mm	Ø 230 x 300D	Ø 230 x 410D	Ø 250 x 450D	Ø 325 x 410D	Ø 325 x 600D
Outer dimensions (WxDxH)	mm	570 x 500 x 450	570 x 610 x 450	570 x 775 x 450	665 x 705 x 450	665 x 895 x 545
Loading type	-	Front Loading	Front Loading	Front Loading	Front Loading	Front Loading
Heater power	kW	2 KW	2.4 KW	2.4 KW	2.8 KW	3.2 KW
Trays	nos	2	2	3	3	3
Special Cycles	-	Vacuuum / Helix / BD Test	Vacuuum / Helix / BD Test	Vacuuum / Helix / BD Test	Vacuuum / Helix / BD Test	Vacuuum / Helix / BD Test
Operating Temp. / Pressure	°C / kPA	121°C / 134°C	121°C / 134°C	121°C / 134°C	121°C / 134°C	121°C / 134°C
Duty Cycle	-	8 cycles / Day	8 cycles / Day	8 cycles / Day	8 cycles / Day	8 cycles / Day
Min Installation Size	mm	600 x 470	650 x 470	790 x 470	725x 570	910 x 570
Total Cycle Time	min	20 min - 60 min	30 min - 75 min	30 min - 75 min	30 min - 75 min	30 min - 75 min
Cycle Record	-	Thermal Printer	Thermal Printer	Thermal Printer	Thermal Printer	Thermal Printer
Packing Dimensions (WxDxH)	mm	620 x 550 x 490	620 x 660 x 490	620 x 825 x 490	700 x 750 x 600	710 x 950 x 600
Net / Gross Weight	kg	45 kg	52 kg	55 kg	65 kg	85 kg

For Details Contact :

MEDISTER

Infection Control Systems Pvt. Ltd.

Office : G-119, Ground Floor, Shagun Arcade Premises, Near Dindoshi Bus Depot, Film City Road, Malad (East), Mumbai - 400 097.

Phone : +91-22-28495501 / 28497725, Email : admin@medister.in;

Website : www.medister.in

Our Sales / Service Network :

Maharashtra / New Delhi / West Bengal / Tamilnadu / Karnataka / Andra Pradesh / Kerala / Madhya Pradesh / Gujarat / Rajasthan / Punjab / Uttar Pradesh / Orrisa / Assam / Himachal Pradesh / Uttarakhand / Goa / Bihar / Jharkhand / Nagpur / Chattisgarh.

Authorised Distributor :